



December 27, 2019

Inmode MD Ltd.  
% Amit Goren  
Regulatory Manager  
A. Stein- Regulatory Affairs Consulting Ltd.  
20 Hata'as Str., Suite 102  
Kfar Saba, 4442520 Il

Re: K192695

Trade/Device Name: InMode System with the Morpheus8 (Fractora) Applicators  
Regulation Number: 21 CFR 878.4400  
Regulation Name: Electrosurgical Cutting And Coagulation Device And Accessories  
Regulatory Class: Class II  
Product Code: GEI  
Dated: November 25, 2019  
Received: November 27, 2019

Dear Amit Goren:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For     Long Chen, Ph.D.  
         Assistant Director  
         DHT4A: Division of General Surgery Devices  
         OHT4: Office of Surgical  
                and Infection Control Devices  
         Office of Product Evaluation and Quality  
         Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration

## Indications for Use

Form Approved: OMB No. 0910-0120  
Expiration Date: 06/30/2020  
See PRA Statement below.

510(k) Number (if known)

Device Name

InMode System with the Morpheus8 (Fractora) Applicators

Indications for Use (Describe)

The InMode System with the Morpheus8 (Fractora) Applicators is intended for use in Dermatological and General Surgical procedures for Electrocoagulation and Homeostasis.

At higher energy levels greater than 62 mJ/pin, use of the Morpheus8 (Fractora) Applicators is limited to Skin Types I-IV.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASaff@fda.hhs.gov](mailto:PRASaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**510(K) SUMMARY**  
**INMODE RF SYSTEM**

**510(k) Number K192695**

**Applicant Name:**

Company Name: InMode Ltd  
Address: Tabor Building, Shaar Yokneam  
POB 44, Yokneam 20692  
Israel  
Tel: +972-4-9097470  
Fax: +972-4-9097471  
E-mail: [amit@asteinrac.com](mailto:amit@asteinrac.com)

**Contact Person:**

Official Correspondent: Amit Goren  
Company Name: A. Stein – Regulatory Affairs Consulting Ltd.  
Address: 20 Hata'as Str., Suite 102  
Kfar Saba 4442520 Israel  
Tel: + 972-9-7670002  
Fax: +972-9-7668534  
E-mail: [amit@asteinrac.com](mailto:amit@asteinrac.com)

**Date Prepared:** November 25, 2019

**Trade Name:** InMode System with the Morpheus8 (Fractora) Applicators

**Classification Name:** CFR Classification section 878.4400; (Product code GEI)

**Classification:** Class II Medical Device

**Predicate Device:**

The subject device is substantially equivalent to the following predicate device:

| Device                                                 | Manufacturer   | 510(k) No. |
|--------------------------------------------------------|----------------|------------|
| InMode System with the<br>Fractora3D/3D-90 Applicators | InMode MD Ltd. | K180189    |

**Device Description:**

The InMode System with the Morpheus8 (Fractora) Applicators is a computerized, programmed, RF technology-based device intended for dermatological applications which requires skin electrocoagulation and hemostasis.

The device platform is basically constituted on the same system platform as FDA cleared for the proposed predicate device; The InMode System with the Fractora3D/3D-90 Applicators (K180189). The device applicators employ fractional RF multi-electrode technology for procedures requiring electrocoagulation and hemostasis. The Morpheus8 (Fractora) Applicators are designed to deliver radiofrequency energy to the skin in a fractional manner, via an array of multielectrode pins.

The Device provides enhanced safety while minimizing possible side effects by monitoring RF parameters. The InMode System with the Morpheus8 (Fractora) Applicators consists of an AC/DC power supply unit, RF generator, controller and user interface including LCD touch screen. The Morpheus8 (Fractora) Applicators are connected to the console via a cable and a foot switch activates the energy delivery to the applicator. The subject device applicator comprises a disposable, single use, fractional RF electrode pins tip head.

The InMode System with the Morpheus8 (Fractora) Applicators comprises the following applicator tip heads:

- Morpheus8 24 Pin Applicator (Fractora3D) tip (FDA Cleared in K180189)
- Morpheus8 40 Pin treatment tip (New tip)
- Morpheus8 12 Pin treatment tip (New tip)
- Morpheus8 T Pin treatment tip (New tip)

**Device Specifications:**

|                                        |                         |
|----------------------------------------|-------------------------|
| Main Line Frequency (nominal)          | 50 - 60 Hz              |
| Input Voltage (nominal)                | 100 - 240 VAC           |
| Electrosurgical Unit dimensions (inch) | 18.2''W x18.2''D x40''H |
| Platform weight (lb.)                  | 70.5                    |
| Applicator weight (lb.)                | 1.26                    |
| RF Max Output Power (Watt)             | 65                      |
| RF Output Frequency (MHz)              | 1± 2%                   |

**Intended Use/Indication for Use:**

The InMode System with the Morpheus8 (Fractora) Applicators is intended for use in Dermatologic and General Surgical procedures for Electrocoagulation and Hemostasis.

At higher energy levels greater than 62 mJ/pin, use of the Morpheus8 (Fractora) Applicator is limited to Skin Types I-IV.

**Performance Standards:**

The InMode System with the Morpheus8 (Fractora) Applicators complies with the following FDA recognized consensus standards:

- AAMI/ANSI 60601-1 (2012), Medical Electrical Equipment - Part 1: General Requirements For Basic Safety And Essential Performance (IEC 60601-1:2005, Mod).
- IEC 60601-1-2 (Edition 3.0, 2007), Medical Electrical Equipment - Part 1-2: General Requirements For Basic Safety And Essential Performance - Collateral Standard: Electromagnetic Compatibility - Requirements And Tests.
- IEC 60601-2-2 Edition 5.0 2009-02, Medical Electrical Equipment - Part 2-2: Particular Requirements For The Basic Safety And Essential Performance Of High Frequency Surgical Equipment And High Frequency Surgical Accessories.

**Non-Clinical (Bench) Performance Data:**

The InMode System with the Morpheus8 (Fractora) Applicators was tested for its applicator tip RF output performance specifications. The device utilizing the Morpheus8 12, 24, 40 pin and T tip heads preset to different RF energy levels, was tested under different resistance loads with the following parameters measured and compared to the results obtained in similar testing conducted on the predicate device: pulse duration, energy per pin, total energy, and pick to pick voltage. The results demonstrate equivalent measurements obtained for the subject and predicate devices.

**Animal Performance Data / Histology Data:**

Not Applicable

**Clinical Performance Data:**

Not Applicable

## Biocompatibility

All of the modified device materials which come in direct contact with the patient skin are biocompatible and identical to the materials used in the predicate device manufacturing.

## Substantial Equivalence:

A comparison table is provided below comparing the intended use and basic technological characteristics of the InMode System with the Morpheus8 (Fractora) Applicators to the intended use and basic technological characteristics of the predicate device.

| Technological Characteristic                   | InMode System with the Morpheus8 (Fractora) Applicators<br>InMode Ltd.<br>(K192695)                                                                                                                                     | InMode System with the Fractora3D/3D-90 Applicators<br>InMode MD Ltd.<br>(K180189) |
|------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| Dimensions                                     | 46cm W x 46cm D x 100cm H<br>[18.2'' W x 18.2'' D x 40'' H]                                                                                                                                                             | 46cm W x 46cm D x 100cm H<br>[18.2'' W x 18.2'' D x 40'' H]                        |
| Weight<br>Platform weight<br>Applicator weight | 32 Kg (70.5 lbs.)<br>0.4 Kg (1.26 lbs.)<br>Tip weight - 0.02 Kg                                                                                                                                                         | 32 Kg (70.5 lbs.)<br>0.4 Kg (1.26 lbs.)<br>Tip weight - 0.02 Kg                    |
| Applicator Tip Heads:                          | <ul style="list-style-type: none"> <li>• Morpheus8 (Fractora) 12 tip head</li> <li>• Morpheus8 (Fractora) 24 tip head</li> <li>• Morpheus8 (Fractora) 40 tip head</li> <li>• Morpheus8 (Fractora) T tip head</li> </ul> | Morpheus8 (Fractora3D) 24 tip head                                                 |
| Number of pins                                 | 12-40 pins                                                                                                                                                                                                              | 24 pins                                                                            |
| Maximal Treatment depth                        | 4.0mm<br>T tip head – 0.5mm (Fixed)                                                                                                                                                                                     | 4.0mm                                                                              |
| Cable Dimensions:                              | 250 cm                                                                                                                                                                                                                  | 250 cm                                                                             |
| Performance                                    | Frequency: 1 MHz                                                                                                                                                                                                        | Frequency: 1 MHz                                                                   |

| <b>Technological Characteristic</b> | <b>InMode System with the Morpheus8 (Fractora) Applicators<br/>InMode Ltd.<br/>(K192695)</b>        | <b>InMode System with the Fractora3D/3D-90 Applicators<br/>InMode MD Ltd.<br/>(K180189)</b>         |
|-------------------------------------|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
|                                     | Maximal RF output power: 65W                                                                        | Maximal RF output power: 65W                                                                        |
| Standards Met                       | IEC 60601-1<br>IEC 60601-1-2<br>ANSI AAMI 60601-2-2 for safety of high frequency surgical equipment | IEC 60601-1<br>IEC 60601-1-2<br>ANSI AAMI 60601-2-2 for safety of high frequency surgical equipment |
| Reprocessing & Sterilization        | For single use.<br>All applicator tip heads are gamma sterilized.                                   | For single use.<br>Non-sterilized. To be sterilized prior to use.                                   |

### **Comparison Discussion:**

The indications for use and technological characteristics of the modified device are substantially equivalent to the indications for use and technological characteristics of the predicate device.

The design and components of the modified device, including the console (with power supply, RF generator, controller and display panel) and the Applicator handle (with cable and connector to console) are identical to the design and components found in the predicate device. The modified device comprises the predicate device Morpheus8 (Fractora) 24 pin tip head and additional three compatible tip heads; the Morpheus8 (Fractora) 12, 40 pin and T tip heads. The performance specifications (including RF frequency, pulse duration, RF current distribution and RF energy per pin) of the modified device are comparable to those of the predicate device. The safety features and compliance with safety standards in the modified device are identical to the safety features and compliance with safety standards found in the predicate device. Patient contact materials are also identical. Any minor differences in the technological characteristics do not raise new safety or effectiveness concerns. Furthermore, the modified device underwent bench performance testing, including software validation testing, electrical and mechanical safety testing according to IEC 60601-1 and electromagnetic compatibility testing according to IEC 60601-1-2 and IEC 60601-2-2. Sterilization methods applied on the device applicator tip heads were validated in accordance with ISO 11137-1.

The performance tests demonstrated that the minor differences in the device design and specifications meet the system requirements and do not raise new safety or effectiveness concerns.

Consequently, it can be concluded that the InMode System with the Morpheus8 (Fractora) Applicators is substantially equivalent to the predicate InMode System with the Fractora3D/3D-90 Applicators, FDA cleared in 510(k) K180189, and therefore, may be legally marketed in the USA.