

User's Manual

Home-care Beauty Device

Photo PLUS EX eye pro

HRF-20-EYE



Table of contents

Safety Precautions	1
When Using Outside of Japan	6
Features	7
Names and Functions of Components	7
Before Use	8
How to Use	13
After Treatment	17
Troubleshooting	18
After-sales Service	21
●Product specifications	22

Thank you for purchasing YA-MAN product.

- Please read this User's Manual carefully for safe and proper use.
- After reading the contents hereof, make sure to store this User's Manual in an accessible location.
- Make sure to include this User's Manual if you are transferring ownership of the device to another person.

Safety Precautions

- Carefully read these Safety Precautions before using the device.
- Use only as directed.
- Immediately discontinue use if any malfunction or unusual reaction to treatment occurs, and ask for a repair.
- Do not use the device in a defective state, as doing so could cause physical harm to the user.
- Precautions listed in the User's Manual are very important. For your safety, please follow them carefully.
- Precautions are classified as DANGER, WARNING, or ATTENTION in order to clearly specify the level of danger and possible injury or damage that may occur if the device is misused.

● Classifications of Precautions

	DANGER	Improper use of this device poses a severe and likely risk of death or serious injury.
	WARNING	Improper use of this device may pose a risk of death or serious injury.
	ATTENTION	Improper use of this device may result in fire, injury, or property damage. * "Property damage" refers to expanded damage affecting a dwelling, household belongings, or domestic animals and pets.

● Symbols

	This symbol shows you what you must not do.
	This symbol shows you what you must do.
	This symbol shows you what you must be careful with.



DANGER

Do not use under any of the following conditions.



Prohibited

Use of the device could cause an electronic medical device to malfunction.

◆ Persons using electronic medical devices.

- Pacemaker or other implanted electronic medical device.
- Artificial heart-lung machine or other electronic life-support medical device.
- Electrocardiograph or other electronic medical devices worn on the body.



WARNING

The following persons should not utilize the device.

May cause an accident or other problems.

- Those undergoing treatment by a doctor, those in a physical growth stage, those with acute medical conditions, tubercular conditions,
- fever, high blood pressure, communicable diseases, malignant tumors, blood diseases;
- Those with diseases of internal organs, skin conditions, allergies, atopic dermatitis,
- Those who are experiencing skin problems or severe fatigue.
- Those with sensitive or delicate skin; in poor health.
- Those who are intoxicated; during their pregnancies, menstruations or breast feeding periods
- Those who take medication that induces drowsiness.
- Those who have unhealthy skin blemishes caused by physical irritation.
- Those who have skin irritation or redness caused by sunburn.
- Those who have capillary dilatation due to long term use of steroid hormonal drugs or a disorder of liver function.
- Those who had an operation on their eye(s) or faces; who have a disease of the eye(s) or face; those who are currently undergoing treatment for a disease of the eye(s) or face.
- Those with excessive myopia (-6D or higher).



Prohibited

Do not use the device on various sites as indicated below and outside of specified sites on the face.

May cause an accident or other problems.

- Body parts with cuts or that have undergone cosmetic surgery; inside the mouth; in the pubic region; on mucous membranes; acnes; areas inflamed due to suppuration; areas with itchiness or inflammation.
 - Swollen moles, viral warts, on the Adam's apple, on the eyeballs, eyelids (except when using Cool Mode)
- Around the eye area (when not using the device with Eye Care or Cool Modes)
* RF LED Mode can be used over eyelids and around the eyes.
When using the LED light on RF LED Mode on the eyelids and areas around the eyes, close your eyes and hold the device 5cm away from the eyes.
Do not look directly into the LED light.

If you are not sure, consult your doctor before use.



Prohibited

Do not apply over the same area continuously.
Always move the device over the skin while treating. Do not use on body parts that cannot be reached by the entire head surface.
May cause a burn injury or a low-temperature burn injury



Prohibited

Do not put in fire; heat; charge, use or leave in a hot place.

The device features a built-in rechargeable battery. May generate heat, ignite, or explode.



User restrictions

Keep away from children and do not allow children to use the device.

Failure to comply may cause an accident or injury.



Non-intended uses are prohibited.

Do not use the device for any purpose other than those stated in the User's Manual.

Failure to comply may cause an injury or a failure of the device.



WARNING



Prohibited

Do not use the device near a flammable substance (such as alcohol, thinner, benzene, nail polish remover, or spray). Failure to comply may cause a fire, electric shock or injury.



Overloading an electric circuit prohibited.

Use only with the supplied AC adapter.
Directly plug the AC adapter into a wall outlet and Do not overload an electric circuit.
Failure to comply may cause an electric shock, short circuit, or fire.



User restrictions

Do not allow persons who are unable to express their intentions on their own to use the device.
Use by unattended persons with physical disabilities should be avoided.
Failure to comply may cause an accident or injury.



Prohibition against use when an odor or abnormal generation of heat is detected

If any problem occurs during use, including abnormal odor or excess heat, turn off the power.
May cause the device to become damaged or to ignite.



Must comply

Check to make sure that there is no water on the adapter plug or power outlet before charging.
Failure to comply may cause an electric shock, short circuit, or fire.



Careless or reckless use is prohibited.

Do not cut, damage, modify, excessively bend, pull, twist, or wrap the AC adapter.

Placing a heavy item on top of the AC adapter or pinching them between other items, could damage them and cause fire or electric shock.



Use of a damaged cord or plug prohibited.

Do not use if the AC adapter is damaged or the wall socket is loose.

Failure to comply may cause an electric shock, short circuit, or fire.
Contact your local distributor for advice on repairing the device.



Do not use near water

Be sure that no water gets on the unit or AC adapter cord. (Do not apply water to or wash the device in water.) Also, do not use near water, in the bathroom, or while bathing.

Failure to comply may cause an electric shock, short circuit, or fire.



Disassembly prohibited

Do not disassemble, attempt to repair, or modify this device in any way.
Failure to comply may cause a fire, electric shock or injury.
Contact your local distributor for advice on repairing the device.



Must comply

Be sure to turn off the power when not in use or when cleaning the unit to avoid a potential electric shock, short circuit, or the fire hazard.
Failure to comply may cause an electric shock, short circuit, or fire.



WARNING



Must comply

Use in an operating voltage range from 100 to 240 volts. Make sure that the AC adapter is fully inserted into an electrical outlet.
May cause a fire, electric shock or accident.



Disconnecting the AC adapter.

Unplug the AC adapter from the power outlet when not in use. Do not insert or disconnect the plug with wet hands.
Failure to comply may cause an electric shock, short circuit, or fire.



Disconnect the plug

When unplugging the AC adapter from the unit or a power outlet, do so by gripping the plug. Do not insert or disconnect the plug with wet hands.
Failure to comply may cause an electric shock, short circuit, or fire.



Combustion warning

Do not leave in hot or humid places, such as near flame, in bright sunlight, or in a bathroom. Do not bring the device close to a heater.
May generate heat, ignite, or explode.



Disconnecting the AC adapter.

In the unlikely event that an abnormality arises, immediately unplug the AC adapter.
Do not use the device if smoke, an abnormal noise, or an odor is detected or if it is otherwise operating in an abnormal state.



ATTENTION



Prohibited

Do not use under any of the following conditions.

Failure to comply may cause an accident, injury, or a failure of the device.

- Do not use the device if your skin is experiencing a rash, inflammation, itchiness, or other abnormalities.
- Do not use if damaged due to dropping or impact with another item.
- Do not use continuously for an overly long period (for example, through commercial use).



Prohibited

Do not look directly into the LED light when using the device on RF LED Mode. When using the LED light to treat areas around the eyes and the eyelids, be sure to close your eyes and keep the electrode head a minimum of 5 cm away from the skin.
Failure to comply may cause injury to your eyes or skin.



CAUTION



Prohibited

Do not wipe the device with alcohol, thinner, benzene, nail polish remover, or other such substance.

Failure to comply may cause damage to the surface of the device.



Prohibited

Do not use the device on pets.

May result in an injury.



Prohibited

Avoid contact between the probe head and metal surfaces.

May result in a failure of the device and injury.



Dropping prohibited

Do not place the device in an unstable location. Do not drop or strike the device.

May cause a failure of the device. An injury may also be caused if the device is dropped on a foot or other exposed body part.



Prohibited

Do not use Cool Mode successively.

If used for successive treatments, the head will not fully cool. To resume using the device, wait at least 30 minutes before proceeding.



Must comply

If you are currently undergoing dental treatment or have teeth fillings, choose a lower intensity level and apply the device on a different area.

Those who are currently undergoing dental treatment or have dental fillings may experience unpleasant sensations in their teeth due to the ion/EMS stimulation. Choose a lower ion/EMS level and apply the device on a different area.



Must comply

Be sure to comply with the following items.

Failure to comply may cause an accident, injury, or a failure of the device.

- Do not use continuously on the same body part.
- Use the device appropriately based on your skin conditions.
- After sunbathing, do not use the device in combination with peeling, skin packs, or other treatments.
- Always turn off the power after use.



Must comply

Stop using if you feel any physical discomfort during use.

Carefully check that your skin is showing no abnormality during a treatment. Do not use the device if a skin abnormality is perceived.



Must comply

Always select Level 1 upon the first use until you become accustomed to use.

Failure to comply may cause an accident or injury.

When Using Outside of Japan



WARNING



Must comply

Be sure to use in an operating voltage range from 100 to 240 volts.

Failure to comply may cause a fire, electric shock, or malfunction.



Be sure to check the plug configuration and the voltage of the power outlet before using.

Failure to comply may cause the generation of heat, a fire, or a malfunction of the device.

- The plug configuration and voltage varies depending on the region and the type of accommodation. Always be sure to check for compatibility before using the device.
- For information about plug adapters, consult a travel agency or a consumer electronics retailer.

Features

This product is a facial beauty device that provides skin care through Ion Cleansing, Ion Moisture, Micro-current, EMS, Vibration, Photo, and Cool functions while warming the skin through RF Deep Heating (hereinafter referred to as "RF"). It can be used to suit your personal skin condition and to achieve your objective by combining the following available treatment Modes.


- RF: Warming the skin^{*1} to prepare the better condition for the treatment.
- Ion Cleansing: absorbs dirt which cannot be removed with normal face washing with the power of ions.
- Ion Moisturizing: helps skin lotion to penetrate into the skin ^{*1} more effectively with the power of ions.
- Micro-current: treats around the eye area with gentle electrical current.
- EMS: Stimulates the muscles with an electrical current.
- Vibration: Supports a skincare with gentle vibration.
- Photo: Treats your skin with LED Photo (red light) treatment.
- Cooling: Tightens skin with cooling effect.

*1 To stratum corneum


●Cotton Stopper (Silver Ring)

*Cotton stopper is attached to electrode head when purchased.

*Use this stopper to apply cotton to the head.



●AC adapter




(*) Design and specifications of the main unit and accessories are subject to change for improvements without prior notes.

Before Use

Before Use


Names and Functions of Components

●Main unit



How to read the panel

●Panel/Buttons



- Press the Power/Mode button for about three second to turn on the device. The treatment will automatically start with Clean Mode.
- Treatment Mode can be changed by pressing the Power/Mode button during operation.
- For Clean/ Moisture Modes, you can select the Ion level. For EMS Up Mode, you can select the EMS level. Each Mode has 5 intensity levels which can be adjusted with the Level/Vibration button.
- For Clean/Moisture Mode, you can press the Level/Vibration button for approx. 3 seconds to turn vibrations ON/OFF.
- Press and hold the Power/Mode button to turn the power off.
- * The device does not vibrate on Eye Care/EMS Up/RF LED/Cool Modes.
- * Intensity level cannot be selected for Eye Care/RF LED/Cool Modes.
- * RF intensity level cannot be selected.


● Mode indicator

Each time you press the Power/Mode button, the treatment mode will be selected in the following order:

Clean→Moisture→Eye Care→EMS Up→RF LED→Cool→Clean (repeat)

Mode	CLEAN	MOISTURE	EYE CARE	EMS UP	RF LED	COOL
Function	RF/ Vibration/ Ion Cleansing	RF/ Vibration/ Ion Moisture	RF/Micro-current	RF/EMS	RF/ LED	Cool
Mode Indicator	LEVEL   CLEAN	LEVEL   MOISTURE	  MOISTURE	LEVEL   EMS UP	 RF LED	 COOL

* When using the device with Clean/Moisture/ EMS Up Mode, the RF Mode indicator repeats to blink on and off.
The device emits one short beep when the (RF) function is operating (RF Mode will be displayed on the Indicator panel).
The device emits two short beeps when the Ion or EMS function is operating (RF Mode will not be displayed on the Indicator panel).



● Level indicator

- Press the Level/Vibration button to select Ion input/output and EMS levels.
- Select from 5 intensity levels. Please select the intensity level which you feel the most comfortable.



- * You can select stronger/weaker ions for Clean and Moisture Modes.
- * You can select stronger/weaker EMS level for EMS Up Mode.
- * Always select Level 1 upon the first use until you are accustomed to use.
- * Intensity level cannot be selected for Eye Care/RF LED/Cool Mode.
- * RF intensity level cannot be selected.


● Charging indicator

The charging indicator will appear on the indicator panel and will flash while the unit is charging.

When charging has completed, the charging indicator will turn solid.

* You can use the device while charging with the AC adapter.

* Charging indicator will not appear during use.



● Standby/Auto-off (Automatic Power Off)

When treatment on Clean/Moisture/Eye Care/EMS Up/RF LED Mode has finished, the device will emit three short beeps. The indicator for the next mode will blink for approx. 2 minutes in Standby Mode (device is not operating).

When you press the Power/Mode button, the next Mode which is waiting in Standby Mode will start.

If the next mode is not started within approx. 2 minutes, the device will emit two short beeps and the power will turn off automatically (Auto-off).

When treatment on Cool Mode has finished, the power will switch off automatically (Auto-off) without the device going to Standby Mode. (The device emits two short beeps and the power switches off.)

● Vibration

- With Clean/Moisture Mode, the vibration function can be turned on/off by pressing the Level/Vibration button for approx. 3 seconds or more.

* Vibration rhythm switches automatically.

* Vibration does not operate with Eye Care/EMS Up/RF LED/Cool Mode.

Pre-testing on skin

* When using the device or a cosmetic product for the first time, make sure to perform a pre-test to check its compatibility with your skin.

- Clean your skin.**
 - Wash your arm with soap and water and wipe dry with a towel.
- Apply cosmetics (gel or lotion) or water to your arm and use the device for approx. three minutes on Moisture Mode.**
 - Do not wet or rub the pre-testing site.

If you experience any skin abnormality (redness or burning sensation) at the pre-testing site within 24 hours, stop to use the device.



- Make sure to perform a pre-test in order to prevent skin problems from occurring.
- If a rash, reddening, itch, stinging sensation, or other abnormality occurs at the pre-testing site, immediately stop the pre-test and wash the pre-testing site.
- During treatment, carefully check that your skin is not showing any abnormality. Stop using the device immediately if the treatment does not agree with your skin.

Charging the Battery

The product is rechargeable and can be used either with or without connecting to an AC adapter.

Always use the device after charging when you use device without an electrical cord.

- Before using for the first time, or if you have not used the device for an extended period, fully charge the battery and make sure that the Charging indicator has turned from flashing to solid, indicating a full charge.

●Charging time

Charge time	Operating time
Approx. 2.5 hours	Approx. 40 minutes

- * The time indicated at left is when the device is charged at an ambient temperature of around 25°C.
- * Charging time and operating time may vary depending on the operating Mode, operating level, and/or operating environment.
- * Actual charging time will differ depending upon the remaining battery and the charging environment.
- * Operating time is how long it can be used immediately after a full charge.

●How to Charge

1 Insert the AC adapter plug into the adapter socket of the main unit.

2 Plug the AC adapter into the wall socket.

- Use only with the supplied AC adapter.
- Directly plug the AC adapter into a wall outlet and do not overload an electric circuit .
- When charging had started, the Charge indicator will flash. When charging is complete, the Charge indicator will turn solid.

3 When charging is complete, unplug the AC adapter from the wall socket and then remove the adapter plug from the main unit.

About the Treatment Modes


Mode	Product Description	Function	Recommendations	Time	Level	Frequency (rough guide)
CLEAN	Absorbs dirt which cannot be removed with normal face washing and tone your skin, while deeply warming the skin ¹ .	· RF · Ion Cleansing · Vibration	RF Booster Pad Pore Refining Deep Cleansing (or cotton + your regular skin lotion ²)	Approx. 3 minutes	Ion: 5 levels	2 to 3 times a week
MOISTURE	Helps beauty ingredients to penetrate into the stratum corneum and moisturize your skin, while deeply warming the skin ¹ .	· RF · Ion Moisture · Vibration	RF Booster Pad Clear Deep Moisturizer (If moisture is insufficient, add skin lotion remaining in the pouch.)	Approx. 3 minutes	Ion: 5 levels	
EYE CARE	A gentle electrical Micro-current treats the area around the eyes to improve elasticity.	· RF · Micro-current		Approx. 2 minutes	N/A	5 times a week
EMS UP	Stimulates the muscles to make your skin elastic while deeply warming the skin ¹ . *Can also be used on the neck and shoulders.	· RF · EMS		Approx. 3 minutes	EMS: 5 levels	2 to 3 times a week
RF LED	Treats the skin with light while deeply warming the skin. ¹	· RF · LED		Approx. 3 minutes	N/A	
COOL	Improves skin elasticity and skin texture with cooling treatment.	· Cooling	Place on the skin directly.	Approx. 3 minutes	N/A	Everyday use is possible

- Be sure to move the unit on your skin while using it. (except when using the device on Eye Care/Cool Mode)
- Use of cottons on modes other than Clean/ Cool Mode is optional.
- * The device may get hot after using the Cool Mode.
Do not use Cool Mode successively.
Allow the device to rest (with power off) at least half an hour after use.
- * 1 To stratum corneum
- * 2 If using your regular skin lotion, be sure that the skin lotion is water soluble.
Do not use oil-based or scrub-type skin lotions, as this could cause the device to malfunction or skin problems.

Attaching Cotton to the Head

● For how to use RF Booster Pads, please refer to the instructions provided separately.

1 Remove the cotton stopper from the electrode head.



2 Place the cotton on top of the electrode head and press the cotton stopper down over the cotton.

- * If the cotton stopper does not fit properly on the electrode head, use a thinner cotton.
- * If the cotton is too thick, the thermal action of RF treatment may be adversely affected and you will feel less heat.

3 Apply your regular skin lotion or gel onto the cotton.

- * Make sure that the lotion or gel coats the entire surface of the treatment head.
- * If using a smooth (low-viscosity) skin lotion, the feeling of warmth may decrease.
Recommended skin lotions: RF Booster Pad/high-viscosity skin lotions
- * This device is not water-proof so be sure to thoroughly wipe off any lotion or gel left on the device with a dry cloth after each use.
Do not wash with water or other liquids.

How to Use (General Instructions)

● Precautions when using

- Be sure to move the unit on your skin while using it. Using the device continuously on the same area may cause burns or other injuries.
- The temperature of RF may feel different depending upon the condition of your skin and types of cosmetics used together.
- Do not use on eyelids or around the eyes (except when using the device on Eye Care or Cool Mode)
 - * RF LED Mode can be used on eyelids and around the eyes. Close your eyes and hold the device 5cm away from your eyes.
- Do not look directly into the LED light when using the device on RF LED Mode.
- Be sure to hold the device so that your hand is touching the electrode head.
- Recommended ambient temperature is 15°C ~ 30°C.

● Treatment method

1 Remove makeup and gently cleanse your face.

- Remove oil and dirt from your skin so that it is in a clean condition prior to treatment.
- * Cool Mode can be used over makeup.

2 Attach cotton if necessary and apply lotion onto the cotton.

3 Press the Power/Mode button for approx. 3 seconds to turn on the power. Select mode and intensity level.

- Select the mode using the Power/Mode button.
- Select the intensity level using the Level/Vibration button.
- Press and hold the Level/Vibration button for about three seconds to switch the vibration function ON/OFF.
(When using Clean/Moisture Mode,)
- * The device does not vibrate on EYE CARE, EMS UP, RF LED or Cool Mode.
- * Intensity level cannot be selected for Eye Care/RF LED/ Cool Mode
- * RF intensity level cannot be selected.
- * Always select Level 1 upon the first use until you are accustomed to use.

4 Hold the device so that you are touching the hand electrode (silver rim) on the rear of the device.

5 Gently apply the treatment head to your skin and move the head over the area being treated. (Refer to p.15-16 for more details)

- When the Mode is finished, the device will go into standby for 2 min. After using the Cool Mode, the device will automatically turn off.
Either press the Power/Mode button to start the next mode, or press the Power/Mode button for approx. 3 seconds to turn off the power.
- * During use, make sure that the electrode head is firmly touching the skin.

6 Press the Power/Mode button for approx. 3 seconds to turn off the power.

- When treatment on Cool Mode has finished, the power will turn off automatically.

7 After treatment, condition your skin as necessary.

- Please apply skin cream as appropriate.

8 Remove cotton from the electrode head and wipe off any residue of skincare lotion or gel remaining on the electrode head or cotton stopper.

- * This product is not water-proof so be sure to thoroughly wipe off any lotion or gel left on the device with a dry cloth after each use.
Do not wash with water or other liquids.

How to Apply Treatment

RF Treatment Method


Recommended speed of moving the unit: 2cm/second

- * If the temperature feels too hot, move the head more quickly.
- * Make sure the entire surface of the treatment head touches the skin during use.
- * Always move the device over the skin during use. Using the device continuously on the same area may cause a burn or a low-temperature burn injury.
- * Do not place the electrode head in contact with eyelids or eye area.

<Clean, Moisture, RF LED Mode>

Move treatment head over the skin in the direction of arrows shown below.

LED light can be shone on the eyelids and areas around the eyes with RF LED Mode. When using the LED light on RF LED Mode on the eyelids and areas around the eyes, remove the RF Booster Pad, close your eyes, and keep the electrode head a minimum of 5 cm away from the skin.



 Areas you should not treat.


 Areas that can only be treated with LED light of RF LED Mode.
* Close your eyes and hold the electrode head 5cm away from your eyes.

<EMS Up Mode>

Move treatment head over the skin in the direction of arrows shown below.

-  Press and move the device in a circular, massaging motion as indicated by the arrows.

 Areas you should not treat.




<EYE CARE Mode>

Press the treatment head repeatedly from the inner corner to the outer corner of the eye as if stamping on the skin (Press→Release→Move treatment head) Do not use the device with Eye Care Mode on the upper eyelid. Press the treatment head down for approx. 3 seconds.

- * The device can be used around the eyes (lower eyelid, outer corner of the eye) with Eye Care Mode.
- * Do not press the electrode head against the skin for longer than 3 seconds.
- * Do not move the electrode head while pressing down against the skin.

The device will emit two short beeps after one minute. After treatment on one eye has finished, continue your treatment on the other eye.



 Areas you should not treat.

It is recommended to apply the device on one side of your face for 2 minutes, a total of 4 minutes on both sides, for those who are concerned about the dry skin around the eyes.

<Cool Mode>

Apply the device as if you are "stamping" onto your skin: press →lift →move (repeat)

- * The device can be used on eyelids and the area around the eyes with Cool Mode.

To keep the product clean, please be sure to clean it after each use.

●How to clean and store the product

- 1 Ensure the power is off and unplug the AC adapter from the main unit.**
- 2 To clean the electrode head, wipe it with a dampened and well-wrung cloth to remove dirt. Wash the cotton stopper with water.**
 - Remove the cotton stopper from the main unit, and wash it with water then wipe dry.
 - Refrain from using an alkali cleaning agent, alcohol, thinner, benzene, nail polish remover, or other such substance.
 - Do not leave the unit with lotion or gel on it.
- 3 Wipe off any dirt from areas other than the electrode head and cotton stopper with a dry cloth.**

Storing

- If you will not be using the device for an extended period, clean the device before storing.
- When the device is not in use, turn off the main body and disconnect the AC adapter from the main body and wall outlet.
- Avoid storing the device near a stove or heater or in a bathroom or other high-temperature/high-humidity location, location where there is exposure to direct sunlight, or unstable location.
- Keep the device out of the reach of children or pets.

Maintenance and inspections

- If you have not used the device for an extended period, inspect it carefully and check that the device works properly and safely.
- If the device is malfunctioning, go to "Checklist for Possible Problems" (page 18 to 19) and attempt to identify the problem before arranging to have the device repaired. If the problem persists, contact your local distributor (see page 21).
- If the device malfunctions, appears damaged, or needs service, please contact your local distributor.

Checklist for Possible Problems

Symptom	Possible Causes	Solutions
Power does not turn on	→ The Power/Mode button is not being held long enough	→ Press the Power/Mode button for approx. three seconds to turn on the device
	→ Insufficient power	→ Please charge the device or connect the main unit to the AC adapter.
Device only operates for a few minutes even after charging	→ Battery life ended (battery can be charged approx. 300 times).	→ Contact local distributor (see page 21). *Do not try to replace the device by yourself as this may cause a malfunction.
	→ You may have mistakenly pressed the wrong button	→ Press the Power/Mode button for about three second to turn on the device again.
	→ The device is in standby or Auto-off status	→ Press the Power/Mode button (see page 10).
Unit does not start charging. (Charge indicator does not flash)	→ The AC adapter is not inserted into the wall socket (connecting too many plugs in one outlet and such)	→ Connect the AC adapter to a wall outlet.
	→ Adapter is not plugged into main unit	→ Plug adapter into the main unit
		→ Reset by pressing Power/Mode and Level/Vibration buttons at the same time, this may start charging.
Main unit is hot	→ Device was used continuously.	→ Turn off and wait at least 1 hour before use

Symptom	Possible Causes	Solutions
Cool Mode is not very cool	Cool Mode was used continuously	Turn off and wait at least 30 min before use
	Device is being used in a temperature environment exceeding the recommended range.	Try using in the environment with the temperature of 15°C to 30°C
Ion/EMS Mode is too strong	Level is set too high	Lower intensity with Level/Vibration button
	Full surface of electrode head is not touching skin	Apply entire surface of electrode head to the skin
Ion/EMS Mode is too weak (you don't feel anything)	Level is set too low	Increase intensity with Level/Vibration button
	The skin lacks moisture or is dry.	Moisturize skin with lotion or gel; moisturize your hand
Cannot feel EMS Mode	Head is not moved along the muscles	Apply treatment head along muscles (p.16)
RF is too hot	Applying electrode head over one area for a long time	Be sure to keep the roller moving over the skin during operation.
	Moving too slowly	Move the device a little more quickly.
	How the heat feels is vary from person to person	Please shorten the usage time.
Cannot feel any heat of RF	The skin lacks moisture or is dry.	Moisturize your skin with lotion/gel first
	You are moving the head too quickly	Move the treatment head more slowly
RF level does not change		You cannot change RF level; multiple level settings are only for Ion and EMS Modes
Cannot feel vibration		Press the Level/ Vibration button for approx. 3 seconds. (Vibration function only operates for Clean/Moisture Modes)

FAQ (Frequently Asked Questions)

Q	A
How much lotion or gel should I use for each treatment?	This device is designed to achieve better effects with sufficient use of lotions and/or gel. Apply lotion or gel on cotton so that the whole surface of the electrode head is wet. Add more lotion or gel when you feel it is not enough.
Can I use this device in the bathroom?	This product is not waterproof. Do not use the device in the bathroom.
Why does my skin blush after use?	This is due to using RF for deep heating. If it bothers you, use the Cool Mode or apply blue ice to cool down your skin.
Why do I feel a tingly sensation on my face and hands?	This is electrical stimulation from Ion Cleansing, Ion Moisture or EMS and is normal. You can adjust to a comfortable intensity level using the Level/Vibration button. Always select Level 1 upon the first use until you become accustomed to use.
Why do I not feel warmth when using Clean/Moisture/EMS Up/RF LED Mode?	Move the electrode head more slowly or press harder against the skin surface. Moisture may be insufficient, so please add skin lotion or gel. RF repeats turning on and off during treatment. When the  symbol appears, RF is on.
I don't feel any sensations when the device is on Eye Care Mode.	I can barely feel any stimulation from the Micro-current. It becomes difficult to feel the warmth when skin lotion is insufficient, when the electrode head has only been pressed against the skin for a short time, or when the cotton is thick.
Why can't I see any dirt on the cotton?	Dirt is removed even though it may not be visible on the cotton. The amount of visible dirt on the cotton can vary between individuals.
The probe head does not warm up when I use RF Mode.	RF is designed to electrically warm the skin tissue, not the probe head itself.

●About the rechargeable battery

- The rechargeable battery is a consumable item.
- If the battery only powers the device for an extremely short period after charging, it may be at the end of its useful life. With typical use, the battery can be charged approximately 300 times.
- If you need a battery replacement, contact your local distributor. The battery is not user-replaceable.

●About Disposal

- Follow your local government rules for proper disposal of this product.
- This product uses a lithium ion battery.



●Contact

Contact your local distributor for any question you may have regarding the warranty that applies to the device, repairs of the device, or other matters requiring clarification.

●Product Specification

Voltage	DC9V 2.0A
Power Consumption	Approx.13W
Battery	Lithium ion battery
Charging time	Approx. 2.5hours *May vary depending on the charging environment.
Size	Main unit(with Cotton Stopper): Approx.W44×D57×H190(mm)
Weight	Main unit(with Cotton Stopper): Approx.220g
Material	Main Unit: ABS resin Panel: polycarbonate/acrylic resin Electrode head: PET resin/copper (chrome plated) Hand electrode: ABS resin (chrome plated)
Country of Origin	Japan
Accessories	Cotton Stopper, AC Adaptor
Recommended environmental temperature	15℃～30℃